

AMENDMENTS TO THE SPECIFICATION

IN THE SPECIFICATION:

Page 5

Please amend the Specification on page 5 beginning at line 5 as follows:

The Mer-11107 strain was deposited as FERM P-18144 at the National Institute of Bioscience and Human-Technology Agency of Industrial Science and Technology (1-3, Higashi 1-chome Tsukuba-shi, Ibaraki-ken 305-8566 Japan), now reorganized to International Deposit FERM BP-7812 at International Patent Organism Depositary (IPOD) National Institute of Advanced Industrial Science and Technology (Tsukuba Central 6,1-1, Higashi 1-Chome, Tsukuba-shi, Ibaraki-ken 305-8566 Japan), as of December 19, 2000, and then transferred to International Deposit FERM BP-7812 at International Patent Organism Depositary (IPOD) National Institute of Advanced Industrial Science and Technology (Tsukuba Central 6, 1-1, Higashi 1-Chome, Tsukuba-shi, Ibaraki-ken 305-8566 Japan) as of November 27, 2001.

Page 25

Please amend the Specification beginning on page 25 line 20 ending on page 26 line 4 as follows:

The 11107B substance as the substrate may be added in a culture broth or a suspension solution of mycelia either as a powder as it is or as a solution dissolved in a water-soluble solvent, for example, ethanol, methanol, acetone or dimethylsulfoxide. The amount of the 11107B substance added is preferably 50 to 5000 mg per 1 L of the culture broth in the case of the culture broth. After the addition of the substrate, procedures such as flask shaking or tank culture are conducted at 20 to ~~31°C~~ 40°C for about 1 to 5 days to run a reaction in an aerobic condition, whereby the 11107B substance as the substrate can be converted into the 11107D substance.

Page 27

Please amend the Specification on page 27 beginning at line 5 as follows:

One loopful of the slant culture (ISP-2 medium) of Streptomyces sp. Mer-11107 strain (FERM BP-7812) was inoculated

into a 500 mL Erlenmeyer flask containing 50 mL of seed medium (2.0% of glucose, 1.0% of soybean meal (ESUSAN-MEAT manufactured by Ajinomoto Co. Ltd.), 0.5% of yeast extract (manufactured by Oriental Yeast Co., Ltd.), 0.25% of sodium chloride, 0.32% of calcium carbonate, pH 6.8 before sterilized), and it was cultured at 28°C for two days to give the first seed culture broth. 0.1 mL of the culture broth was inoculated into a 500 mL Erlenmeyer flask containing 100 mL of the same seed medium and it was cultured at 28°C for one day to give the second seed culture broth. The second seed culture broth (800 mL) thus obtained was inoculated into a 200 L tank containing 100 L of a production medium (5.0% of soluble starch, 0.8% of Pharmamedia, 0.8% of gluten meal, 0.5% of yeast extract and 0.1% of calcium carbonate, pH 6.8 before sterilized) and it was cultured for five days with ~~flowing air and stirring in~~ the following conditions, to give a culture broth.